262-57-HMS

Position the needle at 45 degrees to the point of maximal pain of the lateral epicondyle. After puncture of the skin, angle the needle parallel to the skin and insert it towards the point of maximal pain on the lateral epicondyle. Inject half the contents as the needle is withdrawn to the skin without exiting the skin.

Rotate the needle 180 degrees (opposite direction) and insert the needle parallel to the skin towards the point of maximal pain on the lateral epicondyle. Inject the remaining contents on withdrawing the needle. Remove the needle from the skin.

Flex and extend the elbow five times and then internally and externally rotate five times.

Rotator Cuff Tendinopathy

- one 1.2ml peri-articular injection of SportVis™ into the subacromial space of the shoulder just above the tendon followed by a second peri-articular injection after 14 days, 22 gauge needle is recommended.

Seat the patient in an upright position, arm relaxed at the side and externally rotated. Use of an ultrasound probe positioned on the lateral shoulder directed in the plane of the supraspinatus tendon to guide the injection is recommended.

Locate the acromion, greater tubercle, head of the humerus and subacromial cleft. Introduce the needle into the cleft 1cm posterior and 2cm distal to the antero-lateral acromial edge. Advance the needle horizontally and in a partly medial direction under the acromion process. When no resistance to the plunger is felt inject SportVis™ over the head of the humerus into the subacromial space taking great care not inject into the tendon.

Warnings and Precautions:

SportVis™ should only be injected by a healthcare professional trained in the procedure.

SportVis™ pre-filled syringes are single use. The contents of the syringe should be used for one injection only. Any remaining sodium hyaluronate should be discarded. If a syringe is retained for a subsequent injection there is a risk of contamination resulting in the possible infection of the patient and/or foreign body reaction.

SportVis™ should not be re-sterilised as the device performance may be compromised which could cause serious harm to the patient's health and safety.

SportVis™ must not be injected into blood vessels because sodium hyaluronate has the potential to occlude the vessels, which could result in embolism or infarction.

SportVis™ should not be injected into a haematoma.

Direct injection into tendons should be avoided as this can lead to rupture. Use of an ultrasound probe to guide the injection will minimise this risk.

Do not inject into the soft tissue of patients if the area of the injection site is infected or where there is evidence of acute or chronic skin disease.

Sodium hyaluronate is manufactured by fermentation of *Streptococcus equi* and rigorously purified. However, the healthcare professional should consider the immunological and other potential risks that can be associated with the injection of any biological material.

There is a risk of infection at the injection site as with any peri-articular procedure.

There is no evidence of the safety of SportVis™ in human pregnancy and lactation.

The safety and effectiveness of SportVis™ has not been tested for children below 18.

Follow national or local guidelines for the safe use and disposal of needles. Obtain prompt medical attention if injury occurs.

Do not use if packaging has been damaged.

Do not use after the expiry date.

Adverse Reactions:

Mild erythema that should resolve with time.

Incompatibilities:

SportVis™ has not been tested for compatibility with other substances for peri-articular injection. Therefore the mixing or simultaneous administration with other peri-articular injectables is not recommended.

Storage:

Store between 2°C and 25°C. Do not freeze. Protect from light.

Do not use if sterile packaging has been damaged.

Do not use after expiry date.

Symbols used on the packaging:

LOT	Lot number
	Expiry date
	For single use only
	Store between 2°C and 25°C
	Protect from light
	Do not freeze
	Do not use if sterile packaging has been damaged
STERILE A	SportVis™ pre-filled syringes are aseptically filled with sterile filtered sodium hyaluronate.
	Manufacturer
	Date of Manufacture
	Full instructions, warnings and precautions are included in this instruction for use leaflet.
Issue Date 04/2013	

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